How to use the Brandel harvester:

- 01. Soak the filter in 0.3% PEI for 20 min
- 02. Make sure that the valve on the waste outlet is closed
- 03. Place your filter and lock the harvester
- 04. Connect the wash water line
- 05. Open the orange cap on the water tubing near the pump
- 06. Turn the harvester on (little switch located at the back left)
- 07. Turn on the vacuum pump (the switch is at the back right)
- 08. Wait a few minutes for vacuum to build
- 09. Place your block
- 10. Switch lever to "Harvest"
- 11. Bring the needles about an eighth of an inch into the block
- 12. Press the black button at the front of the harvester until the water fills all the way to the needles, and you start hearing water being aspirated.
- 13. Bring the block up towards the needles until all liquid is aspirated
- 14. Bring the block back down and repeat two more times.
- 15. Let the water flow until there's no more flow. Wait until vacuum is stable
- 16. Switch lever to "Purge" and let all the water drain. Wait until vacuum is stable.
- 17. Switch lever to "Off"
- 18. Remove your filter and bring it to your bench together with the block
- 19. Take a wet paper towel and place it in the harvester where the filter usually is. Make sure to lock the harvester.
- 20. Switch to the washing block and wash three times. The last wash should be a long one, ~30 sec
- 21. Disconnect everything and turn everything off in the reverse order
- 22. Open the waste outlet to drain the waste into the sink